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### GUARANTEE DOCUMENT

Thank you for purchasing this COMELEC product, this device has been checked in all its components meticulously, so we guarantee the correct operation, but for this to be the case we need that YOU read the instruction book carefully, as this will contribute to the best performance of this article.

P.G.E., S.L guarantees to the buyer of this product, a total guarantee of TWO YEARS, governed by the established conditions whose regulation is found in Royal Legislative Decree 1/2007 of 16 December, which approves the revised text of the General Law for the Defense of Consumers and Users, (articles 114 to 127)

P.G.E., S.L is committed to repairs are completely free, labor and spare parts, including.

The Technical Assistance Service will resolve any repairs that may be required by this apparatus.

Carefully review the INSTRUCTION BOOK included, on the use, maintenance and conservation of the device.

Any anomaly that could occur due to the non-observance of such instructions will not be protected for the guarantee.

The date of sale of the device, for the purposes of warranty period, will be that which appears on the purchase invoice.

Guarantee operation: As regards all those compliance defects detected by the technical assistance, which are manifested within six months of the delivery of the goods, it is assumed that they already existed on that date, the intervention will be completely free, (unless this hypothesis is incompatible with the nature or nature of the conformity defect).

During the following eighteen months, it must be the user who demonstrates. That the lack of conformity it was coming at the time of delivery.

MANUFACTURED CHINA FOR COMELEC

### DECLARATION OF CONFORMITY

Manufacturer Name. PRAGMATICA GRUPO EUROPEO, S.L  
Address: Tapissers 10 Poligono Industrial Els Mollons  
Post Code, 46970 City: Alaquas VALENCIA

We declare under our responsibility that the product:

Product Description: Infrared thermometer

Model : TN 400

Brand: COMELEC

It complies with the European Union legislation in compliance with the 2004/108 / EU guidelines

The standards apply:

EN 55014-1: 2006+A1:+A2:

EN 55014-2: 1997+A1:+A2:

EN 61000-3-2:2006+A1:+A2:

EN 61000-3-3:2013

Name: Pablo Gonzalez

Position : Administrador

Date: 10 Enero 2.019

Signature:

### INFORMATION ON RECYCLING

Disposal of the old appliance, the European Directive 2002/96 / EC on waste equipment Electrical and Electronic (WEEE), establishes that old appliances can not be thrown as municipal waste without normal classification. Appliances have to be collected separately with the aim to optimize the recovery and recycling of the materials they contain and reduce the impact on human health and the environment.

The symbol "container with wheels" crossed out on the product reminds you of obligation, that when you dispose of the device must be collected separately. Consumers should contact the local or retail authority to obtain information about the correct disposal of your old device.

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DATE OF SALE

Seller's signature and stamp

IMPORTANT: To make the guarantee effective, the SELLER must fill out, sign and stamp this Certificate. To receive assistance under the conditions set forth in this certificate, please submit to the Technical Service Authorized, this document next to the Purchase Invoice.

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


2. Buttons: Figure 5

- (1) When pull the trigger , LCD display reading with SCAN icon. Release the trigger , display reading with HOLD icon for 12 sec (approx). Built-in 12 sec auto power off function.
- (2) Laser on/off button
- (3) Celsius / Fahrenheit switch button
- (4) Backlight will be powered off without any operation after 12 seconds

G. Maintenance

- 1. Lens Cleaning:  
Blow off loose particles using clean compressed air. Gently brush remaining debris away with a moist cotton swab. The swab may be moistened with water .
  - 2. Case cleaning: Clean the case with a damp sponge/cloth and mild soap.
- Note:
- 1) Do not use solvent to clean plastic lens.
  - 2) Do not submerge the unit in water.

H. SPECIFICATION

|                        |                                                                                                                                                    |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| TEMPERATURE RANGE      | -50°C ~ 400°C (-58°F ~ 752°F)<br>-50°C ~ 600°C (-58°F ~ 1112°F)                                                                                    |
| ACCURACY               | 0°C ~ 400° C (32°C ~ 752°F) 0°C~ 600°C (32°F~1112°F)<br>± 1.5°C (±2.7 °F) or ±1.5%<br>-50 °C~0°C (-58°F~32°F) ±3°C (±5°F)<br>Which ever is greater |
| RESOLUTION             | 0.1° C O 0,1° F                                                                                                                                    |
| REPEATABILITY          | 1% reading o 1°C                                                                                                                                   |
| RESPONSE TIME          | 500 m Sec, 95% response                                                                                                                            |
| SPECTRAL RESPONSE      | 8-14 um                                                                                                                                            |
| EMISSIVITY             | 0.95 Preset                                                                                                                                        |
| DISTANCE TO SPOT SIZE  | 12: 1                                                                                                                                              |
| OPERATING T EMPERATURE | 0 ~ 40 °C (32 ~ 104°F)                                                                                                                             |
| OPERATING HUMIDITY     | 10 ~ 90% HR non-condensing up to 30°C (86°F)                                                                                                       |
| STORAGE T EMPERATURE   | -20 ~ 60 °C(-4 ~ 140 °F)                                                                                                                           |
| POWER                  | 1.5V AAA * 2 batería                                                                                                                               |
| TYPICAL BATTERY LIFE   | Laser o f f : 12 hrs                                                                                                                               |

A. INTRODUCTION

This infrared thermometer is used for measuring the temperature of the object's surface, which is applicable for various hot, hazardous or hard to reach objects without contact safely and quickly. This unit consist of Optics, Temperature Sensor Signal amplifier , Processing circuit and LCD Display. The Optics collected the infrared energy emitted by object and focus onto the Sensor . Then the sensor translates the energy into an electricity signal. This signal will be turned out to be digital shown on the LCD after the signal amplifier and processing circuit.

B. WARNING & C AUTIONS

1. Warning:

To avoid the potential situation may cause harm or damage to people, please pay attention to the following items:

- 1) Before you use this unit, check on the plastic housing carefully . If there is any damage, do not use it.
- 2) Do not point laser directly at eye or indirectly off reflective surfaces.
- 3) Do not use this unit in the environment of explosive gas, steam or dusty .

2. CAUTION

To avoid the damage of the unit or the target, please protect from the following situations:

- 1). EMF (electro-magnetic fields) from arc welders, induction heaters.
- 2). Thermal shock (caused by large or abrupt ambi- ent temperature changes-allow 30 minutes for unit to stabilize before use.


- 3). Do not leave the unit on or near objects of high temperature

C. Distance to Spot size

- 1. When take measurement, pay attention to the Di- stance to Spot Size. As the Distance (D) from the target surface increases, the spot size (S) of the area measured by the unit becomes larger .

The Distance to Spot size of the unit is 12:1

\*\*\*This unit is equipped with a laser , which is used for aiming.

2. Field of view:

Make sure the target is larger than the unit's spot size. The smaller the target the closer measure distance. When accuracy is critical, make sure the target is at least twice as large as the spot size.

D. Emissivity

Most organic materials and painted or oxidized surfaces have an emissivity of 0.95 (pre-set in the unit). Inaccurate readings will result from measuring shiny or polished metal surfaces. To compensate, cover the target surface with masking tape or flat black paint. Measure the tape or painted surface when the tape or painted reach the same temperature as the material underneath.



E. Operation

1. Operating the unit

- 1). Open the battery door and insert 2\*1.5V AAA batteries properly
- 2). Pull the trigger to turn on the unit;
- 3). Aim at the target surface and pull the trigger , then temperature will be shown on the LCD. This unit is equipped with a laser , which is only used for aiming.


- 2. Locating a Hot Spoto find a hot spot, aim the thermometer outsideof interest, then scan across with an up and down motion until you locate the hot spot.

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
## INSTRUCTIONS MANUAL INFRARED THERMOMETER FOR OBJECT'S SURFACE MODEL TN400

**TECHNICAL SPECIFICATIONS**  
Model : TN 400  
The rated power: 2XAAA 1,5V =  
Power : 1,5mW  
Operation Distance : 3 ~ 5 Mts

**PLEASE READ THE INSTRUCTIONS CAREFULLY OF THIS MANUAL BEFORE STARTING.**

**DO NOT HANDLE THIS INFRARED THERMOMETER INTERNALLY.**

The texts, photos, colors, figures and data correspond to the current technical level from the moment they have been printed. We reserve the right to modifications motivated by the permanent development of the technique in our

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